Non-instru for Use In-sterile Metal Bone Plate

Metal Bone Plate is mainly used for internal fixation of fractures. The locking bone plate is designed and manufactured according to the anatomical characteristics of human bones...

Material: Metal Bone Plate: unalloyed titanium, titanium alloy (Ti6Al4V)
Intended Purpose: Metal Bone Plate is intended to be used with bone screws for temporary fixation, correction and stabilization of limb fractures...

- Indications: 1. Distal humeral fracture. 2. Distal humeral locking plate. 3. Proximal humeral locking plate. 4. Distal femoral locking plate. 5. Proximal femoral locking plate. 6. Distal femoral locking plate. 7. Proximal femoral locking plate. 8. Distal femoral locking plate. 9. Proximal femoral locking plate. 10. Distal femoral locking plate.

- Contraindications: 1. General local infection, osteomyelitis. 2. Severe liver and kidney dysfunction. 3. Mental disease, systemic connective tissue disease. 4. Pregnancy, metastasis or spread of cancer. 5. Fracture of the tibia or fibula. 6. Fracture of the radius or ulna. 7. Fracture of the humerus. 8. Fracture of the femur. 9. Fracture of the tibia or fibula. 10. Fracture of the radius or ulna.

Preparation for use: 1. Cleaning: Use a soft brush to clean the surface of the instrument. 2. Disinfection: Immerse the instrument in a disinfectant solution for 30 minutes. 3. Rinsing: Rinse the instrument with sterile water. 4. Drying: Dry the instrument with sterile gauze.

Storage: 1. Store the instrument in its original packaging. 2. Store in a clean, dry, and well-ventilated area. 3. Avoid exposure to moisture and direct sunlight. 4. Avoid contact with corrosive substances.

Quality Control: 1. Raw material inspection. 2. In-process inspection. 3. Final inspection. 4. Post-market surveillance.

Warranty: 1. The manufacturer warrants the product against defects in materials and workmanship. 2. The warranty is void if the product is used for purposes other than those intended.

Disclaimer: 1. The manufacturer is not responsible for any damage or injury caused by the use of the product. 2. The user should read the instructions carefully before using the product.

Notes: 1. The product is not intended for use in the presence of an active infection. 2. The product is not intended for use in the presence of a tumor.

References: 1. Orthopaedic Surgery, 10th Edition, JB Lippincott Williams & Wilkins. 2. Principles and Practice of Orthopaedics, 10th Edition, JB Lippincott Williams & Wilkins.

Revision History: 1. Revision 1.0: Initial release. 2. Revision 1.1: Updated technical specifications. 3. Revision 1.2: Updated safety information.

Appendix A: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix B: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix C: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix D: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix E: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix F: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix G: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix H: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix I: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix J: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix K: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix L: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix M: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix N: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix O: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix P: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix Q: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix R: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix S: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix T: 1. Diagram of the product. 2. List of components. 3. List of accessories.

Appendix U: 1. Diagram of the product. 2. List of components. 3. List of accessories.

produce negligible artifacts by X-ray, CT, MRI, and ultrasound.
1. Preoperative Plan: (1) Take a X-ray picture of the fractured bone, including adjacent joints. (2) The recovery of the original bone status and the correction of bone malformation are the main purposes of the operation.

2. Surgical Approach: (1) Carry out routine disinfection after successful anesthesia. (2) After the incision, expose to the fracture and lay layer by layer for the functional reduction. (3) The locking bone plate can be followed to implanted after the functional reduction. (4) The locking bone plate is temporarily fixed to the fracture end with bone forceps or kirshner wire. (5) The length of distal and proximal ends of the fracture were measured by a drill depth sounder guided by a guide. (6) Select an appropriate length of the locking screw and screw it into the bone plate with a screwdriver. (7) Complete the installation of the other screws one by one by following the procedure of drilling, depth sounding and locking in locking screw.

(8) Check again if all the screws are tightened. (9) After the fracture fixation is completed, checked again the locking bone plate to see whether the position is appropriate and whether the fracture end achieves the desired stability effect. (10) The whole process should follow the principle of aseptic operation. 3. Removal of Implant: (1) Expose the plate and screws from the original incision and strip tissues on the surface of the bone plate and screws. Use a hexagonal socket screw driver to unscrew. Remove the plate. (2) Disposal of Implant: Implants taken out from body should be handled in an anti-pollution way according to hospital protocol so as to prevent cross infection. Devices must be disposed of as a healthcare medical waste in accordance with hospital procedures. 5. More detailed operation instructions and selection of supporting instruments could be found in the "Operation Manual".

Table with 2 columns: Manufacturer (Double Medical Technology Inc.), Legat Address (No. 18, Shanbenhong East Road, Hainan District, 361026, Xiamen, Fujian, PEOPLE'S REPUBLIC OF CHINA), Telephone Number (+86 592 6087101), SSQP (Summary of safety and clinical performance), Eudamed (https://ec.europa.eu/eudamed/), IVDU (https://www.doublemedical.com)

Note: If any serious incident that has occurred in relation to the device, please report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Notified Body: TÜV BGD Product Service GmbH, Address: Riederstraße 85, 80339 Munich, Germany

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Table with 2 columns: Type of sterilization, Recommended vacuum. Rows include: Minimum sterilization time, Temperature sterilization, Minimum sterilization time, Minimum sterilization time.

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Table with 2 columns: Number of catalog, Manufacturer. Rows include: Date of production, Lot, Authorized representative in the European community, Shelf life, Conforms with instructions, Medical device.

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Organismo Notificado: TÜV SÜD Product Services GmbH Dirección: Riederstraße 85, 80339 Munich, Alemania
Shanghai International Holding Corp., GmbH (Europa) Elberfelder 80, Hamburgo, Alemania