Instructions for Use

Sterile Metal Cannulated Screw

Metal Cannulated Screw consists of a series of screws with different sizes and structures, which is designed according to the anatomical characteristics of human bones. In clinical practice, Metal Cannulated Screw is used alone, which acts as a temporary internal support, provides a stable local environment for the fractured end, and creates conditions for the healing of the fractured end.

I. Material

Metal Cannulated Screws are made of titanium alloy.

II.Intended Purpose

Metal Cannulated Screw is intended for temporary fixation, correction and stabilization of limb fractures.

III.Indication

Metal Cannulated Screw is indicated for internal fixation of limb fractures.

IV. Contraindications

- General or local infection, osteomyelitis.
- Patient can't tolerate operation or anesthesia because of poor health.
- Mental disease, Systemic neurological disease
- Pregnancy, metabolism disorder of calcified tissue
- Incision of the tissue can lead to infection or skin necrosis due to poor local soft tissue conditions.
- Severe osteoporosis prevents effective fixation
- Patients with influence factors for fracture healing (e.g., diabetes, infection, etc.)
- Patients with high degree of comminution, III degree open, severe bone defect and severe damage to surrounding soft tissue.

V.Intended user

The operation should be performed by professional trained surgeons according to the IFU and surgical technique guide strictly. The surgeon should have a thorough understanding of the intended purpose of the product and the surgical techniques used, and should be trained and certified by relevant institutions (e.g., Association for Internal Fixation Research, AO/ASIF).

VI.Patient group
Adults who have the surgical indications and who have no surgical contraindications and can tolerate surgery.

VI.Clinical benefits

Metal Cannulated Screw is used for internal fixation of bone fracture to achieve bone union with expected union rate no less than 90.8%.

VIII. Potential Side Effects

- Implant displacement, screw penetration, loosening or implant fracture
- Nonunion, delayed union or malunion or bone fracture
- Adverse tissue reaction, allergy/hypersensitivity reaction
- Infection
- Soft tissue damage or neuro-vascular damage
- Inflammatory reactions, osteolysis or bone necrosis
- Pain or discomfort
- Impingement or poor joint mechanics
- Wear or corrosion of metal implants.

IX. Warnings and Precautions

1. Responsibility for proper selection of patients, adequate training, experience in the choice and placement of implants and the decision to leave or remove implants postoperatively, rests with the surgeon.

2. The surgeon should discuss the expectations of surgery inherent in the use of the product with the patient preoperatively. Particular attention should be paid to a discussion postoperatively and the necessity for periodic medical follow-up.

3. Suitable implant should be chosen for specific fracture and implant indication. Choose a suitable metal locking bone plate according to preoperative X-ray film. Select a supporting instrument set and metal locking bone screw according to the chosen locking bone plate. Surgery should be finished with the guidance of Manual of Internal Fixation: Techniques Recommended by the AO-ASIF Group and other relevant information.

4. The correct selection of the product is extremely important. The product should be used in the correct anatomic location, consistent with accepted standards for internal fixation. Failure to use the appropriate product for the application may result in a premature clinical failure. Failure to use the proper component to maintain adequate blood supply and provide rigid fixation may result in loosening, bending, or fracturing of the product and/or bone.

5. Careful preoperative planning is a requirement. The surgeons should choose the proper implants according to preoperative images and by using relating measuring instruments such as depth gauges and trials. Instruments which are intended for the implants should be applied in order to obtain a safe combination.

6. In order to avoid or minimize specific risks associated with implantation, surgeons should have been trained, and the operation should be performed strictly according to the surgical technique guide.

7. Implant is not allowed to be reused. The risk of reuse is as follows:

- Metal strength decreases and fracture occurs easily after being implanted.
- As contaminants and bacteria remaining on the surface, re-implantation would result in infection and rejection.
- It can damage the surface anodic oxide film and affect the biological properties of the product.
- Careful handling and storage of the product is required. Scratching or damage to the component can significantly reduce the strength and fatigue resistance of the product.
- Generally, there is no limitation on substance used in the clinical setting, to which the implant might be exposed.
- Once applied, the product should never be reused. Although it may appear undamaged, previous stresses may have created imperfections that could reduce its service life.
- Extreme angles together with small bending radii must be avoided because of the potential risk for postoperative breakage.
- Excessively aggressive use of bending instruments can cause recognizable macroscopic damage to the implant (indentations, elongated screw holes, etc). In such cases, the implant must be exchanged for a new one.
- Unless specifically indicated in the surgical technique guide, the modification of size, shape or surface condition is prohibited after supply.
- Metal Cannulated Screw has not been evaluated for safety in the MR environment. It has not been tested for healing or unwanted movement in the MR environment. The safety of Metal Cannulated Screw in the MR environment is unknown. Performing an MR exam on a person who has this medical device may result in injury or device malfunction.
- Surgeons should verify whether the implants are positioned correctly using adequate imaging technology.

16. Patients should be advised to seek medical opinion before entering potentially adverse environments that could affect the performance of the implant.

17. Premature physical activity or load bearing after surgery is prohibited.

18. The patient should be advised to report any unusual changes of the operated site to his surgeon. The patient should be closely monitored if a change at the fixation site has been detected. The surgeon should evaluate the possibility of subsequent clinical failure, and discuss with the patient the need for any measures deemed necessary to aid healing.

19. Postponed or long-term use may bring about implant looseness or crack; therefore, after bone healing, it's suggested to take out the implant if there is implant crack risk, which may bring medical accident.

20. Considerations for removal of the implant after healing. If the device is not removed after the completion of its intended use, any of the following complications may occur:(1) corrosion, with localized tissue reaction or pain; (2) migration of implant position resulting in injury. (3) Risk of additional injury from postoperative trauma; (4) bending, loosening, and/ or breakage, which could make removal impractical or difficult; (5) pain, discomfort, or abnormal sensations due to the presence of the device; (6) possible increased risk of infection; and (7) bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid refracture. If the patient is older and has a low activity level, the surgeon may choose not to remove the implant thus eliminating the risks involved with a second surgery.

X. Sterilization

The implant is provided as sterile. It is sterilized with γ ray, and the validity period is 5 years.

XI. Labels

	Catalogue number		Manufacturer
	Date of Manufacture		Use -by date
	Batch code		Do not re-use
	Sterilized using irradiation		Do not use if package is damaged
	Authorized Representative in the European Community		Caution
	Keep away from sunlight		Keep dry
	Consult instructions for use		Importer
	Medical device		Double sterile barrier system

XII. Storage
The sterilized product should be stored in a clean environment, protected from direct sunlight and pests.

XIII. Transportation

Avoid collision and compression during storage and transportation.

XIV.Use of Original Products

Implants and instruments are designed to be used together. The use of products from other manufacturers along with Double Medical products can involve incalculable risks, injury or corrosion of the material and misalignment of implant and instruments, impeding functionality, thereby endangering the patient, user or third parties.

XV. Handling Information

Metal Cannulated Screws is made of titanium alloy. The material is biocompatible as widely used in the industry, corrosion-resistant and non-toxic in the biological environment.

XVI.Surgical Technique

1. Preoperative Plan

- Take an X-ray picture of the fractured bone, including adjacent joints.
- The recovery of the original bone status and the correction of bone malformation are the main purposes of the treatment.
- Assess the condition of soft tissue and exam the neurological and vascular functions.
- Treatment of calcaneal osteotomy with 7.0mm HCS

Surgical Approach

(1) patient position: The patient is in a semi-lateral position (Healthy limbs on the bottom, sick limbs on the top)

(2) Surgical approach: Make an incision at the end of peroneal tendon perpendicular to the calcaneal body. Draw the skin apart on both sides with a hook. Care must be taken to protect the peroneal tendon and sural nerve during the incision.

(3) Establishment of osteotomy: Osteotomy incision perpendicular to the calcaneal body with a swing saw should be performed, and medial cortex should be cut with a bone chisel (osteotomy) in order to prevent damage to medial nerve and vascular structure.

(4) Place a guide needle: Transposition the calcaneal body to medial side and insert a guide needle(REF: 112480600, specification:φ2.4×250). Fix the osteotomy end stably in the best position and confirm the position of the guide pin under fluoroscopy.

(5) Measuring depth and selecting screw: Use depth gauge to measure the insertion depth of guide needle, and select cannulated screw of suitable size according to the measuring depth. (REF of depth gauge: 112480400)

(6) Drilling: Place soft tissue protection drill sleeve along the guide pin, use corresponding counter-sink to drill holes in the proximal cortex, and then use corresponding cannulated drill bit to drill holes. The drilling depth is the same as the length of the selected screw.(countersink, type III, REF: 112481900, specification:φ7.5×180, cannulated drill bit, type III, REF: 112481800, specification: 112481800, REF of protection drill sleeve: 112480700).

(7) Screw in: Use cannulated hexagonal screwdriver to screw in the proper position along the guide pin. Confirm the screw position under fluoroscopy, exit the guide pin. Use solid hexagonal screwdriver to finally set for to ensure that the screw tail and bone surface are even (Screw in after other screws in turn according to this operation step).

(8) Suture and dress up the wound.

3. Removal of the implant

Expose the screw from the original incision and remove the surface tissue of the screw. Use a hexagonal screwdriver to remove all the screws. Connect the hexagonal screwdriver and the screw completely so that they are in a straight line in case the tool slips.

Notice:

The handling of gear slip or screw crack

(1) When gear slip occurs and it is impossible for the hexagonal socket screw driver to remove

the screw, connect a slip extractor to an electrical bone drill and insert vertically into the hexagonal socket and the screw can be removed anticlockwise.

(2) When screw crack occurs, connect a crack extractor to an electrical bone drill and drill the bone cortex around the screw anticlockwise and the screw can be removed.

4. Disposal of Implant

Implants taken out from body should be handled in an anti-pollution way according to hospital protocol so as to prevent cross infection.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

5.More detailed operation instructions and selection of supporting instruments could be found in the 《Operation Manual》.

XVII.Manufacturer/ Contact Information

Manufacturer:	Double Medical Technology Inc.
Legal Address:	No. 18, Shanbianhong East Road, Haicang District, 361026, Xiamen, Fujian, PEOPLE'S REPUBLIC OF CHINA
Telephone Number:	+86 592 6087101
SSCP(Summary of safety and clinical performance)	Eudamed(https://ec.europa.eu/tools/eudamed/#/screen/home)
e-IFU	https://www.doublemedicalgp.com

Note: If any serious incident that has occurred in relation to the device, please report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Notified Body: TÜV SÜD Product Service GmbH
Address: Riderstraße 65, 80339 Munich, Germany

	Shanghai International Holding Corp. GmbH(Europe) Eiffestrasse 80, Hamburg, Germany
---	--

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Organismo notificado: TÜV SÜD Product Service GmbH
Shanghai International Holding Corp. GmbH (Europa) Eiffestraße 80, Hamburgo, Alemania