



Anty F-series Femoral Revision Stem System





■ Specification

Diameter	Length				Neck Length	Neck Height	Proximal Outer Diameter				Middle Outer Diameter	Distal Outer Diameter	
	170	205	245	285			170	205	245	285			
1#	14mm	●	●	●	○	42	22	13	14	15	16	14	9
2#	15mm	●	●	●	○	42	22	14	15	16	17	15	10
3#	16mm	●	●	●	○	42	22	15	16	17	17	16	11
4#	17mm	●	●	●	○	42	22	16	17	17	17	17	12
5#	18mm	●	●	●	○	44	24	17	18	17	17	18	13
6#	19mm	●	●	●	○	44	24	18	19	20	21	19	14
7#	20mm	●	●	●	○	44	24	19	20	21	21	20	15
8#	21mm	●	●	●	○	44	24	20	21	21	21	21	16
9#	22mm	×	●	●	○	46	26	×	22	21	21	22	17
10#	23mm	×	○	○	○	46	26	×	23	24	25	23	18
12#	25mm	×	×	×	○	46	26	×	×	×	25	25	20

DOUBLE MEDICAL TECHNOLOGY INC.

No. 18, Shanbianhong East Road, Haicang District, 361026, Xiamen, Fujian, PEOPLE'S REPUBLIC OF CHINA
 Tel: +86 592 6087101 Fax: +86 592 6587078 E-mail: info@double-medical.com Website: www.doublemedicalgp.com

This document is intended for internal reference and use only and is for internal training purposes only. For prohibited contents or precautions, please refer to the instructions. The information contained in this document is not intended in any way to solicit product investment, publicity or promotion, or to use any product beyond its approved indications. It is prohibited to disseminate, copy, publish, copy or quote the content of this document without the express prior written permission of the Company, and it is also prohibited to use any content of this material for any other purpose.